

The HI 839800 thermo-reactor is constructed of durable materials with a vial capacity to perform up to 25 digestions simultaneously. The reactor is equipped with 6 predefined and 2 user-defined temperature profiles. Dedicated temperature profiles for COD and total phosphorous reactions at preset to 150°C. The total nitrogen analysis temperature is pre-set to 105°C. In addition, HI 839800 features LED indicators for visual indication of power on and heating in progress. The heater block temperature can be verified by placing. For more details, please [Contact Us](index.php?option=com_contact&view=contact&id=1) //